

# Certificate of Analysis

## Spectrophotometric Chlorophyll Standards

14 May 2018

Lot CAS158-02

### Original Chlorophyll Solution A:

Mixed: 17 March 2018

Spectrophotometric Data: 9 May 2018

Abs<sub>664</sub> = 0.82646      Abs<sub>750</sub> = 5.2567 E-02

Chlorophyll concentration determination (Jeffrey Method, 1997)

Chlorophyll-a = 8.83 mg/L

Chlorophyll dilution A is provided with this note. The concentration value is guaranteed through May, 2019, but degradation depends on the handling of the standard. Standards must be kept at temperatures below -20° C and in the dark until use. Spectrophotometers should be standardized in subdued light with this solution warmed to room temperature, about 22° C. The determination of chlorophyll-a in aqueous acetone solvent is a sensitive technique. Please pay special attention to the configuration of the spectrophotometer itself. A bandpass of 2.0 nm or less is necessary to properly use this standard.



14 May 2018

John Salinas  
Analytical Chemist

Date



**Reliable Instruments for an Unreliable World**  
1995 N. 1st Street Telephone 408-749-0994  
San Jose, CA 95112 Toll Free 877-316-8049  
www.turnerdesigns.com Fax 408-749-0998