

Trilogy Power Supply Kit P/N 7200-941

This kit is included as an accessory with the Trilogy. It contains interchangeable clips for US, UK, European and Australian wall sockets and voltage stabilization for 10-12V power sources.



Trilogy USB Data Cable P/N 021-7202

This cable is included as an accessory with the Trilogy Model 7200-002. It is a USB cable used for data transfer from the Trilogy to a PC.



Trilogy RS-232 Data Cable P/N 021-0700

This cable is available for users with Model 7200-000. **USB Data Cable P/N 021-7202 will not work.** You must use RS-232 Cable P/N 021-0700. It is a DB9 male to female cable used for data transfer from the Trilogy to a PC.



Communication and Integration Cable P/N 2300-115

This cable can be used as a USB to Serial Adaptor Cable with RS-232 Cable P/N 021-0700. It has been tested using our instruments.

