

Trilogy Power Supply Kit (P/N 7200-941)

This kit is included as an accessory with the Trilogy. It contains interchangeable clips for US, UK, European and Australian wall sockets and voltage stabilization for 10-12V power sources.



Trilogy Data Cable (P/N 021-0700)

This cable is included as an accessory with the Trilogy. It is a DB9 male to female cable used for data transfer from the Trilogy to a PC.



Communication and Integration Cable (P/N 2300-115)

This cable is included as an accessory and can be used as a standard USB to Serial Adaptor Cable. It has been tested using our instruments.

